



Sterilization Container Systems

Surgical Trays & Wire Baskets



Product Catalog

Edition - 01



TABLE OF CONTENTS

06-09

Introduction

10-31

Sterilization Containers

10-20

Standard Models

22-27

Small Models

28-31

Specific Containers

32-51

Perforated Trays & Wire Baskets

32-37

Perforated Trays

38-42

Hybrid Baskets

44-51

Wire Baskets

52-61

Components

63-64

Index



Type and Size Variations

The Sterilization Containers are manufactured in a range of sizes, offering specific solutions to meet the requirements of healthcare professionals.

Designed with user-friendly features that improve convenience in packaging and storage processes, optimizing workflow efficiency for CSSD.



Standard Height Sizes

———— (mm)

100

130

150


200

260

There are 5 different height options available in standard container types

Specific Needs Particular Solutions

The Sterilization Containers are typically available in three standard sizes: full (1/1), medium (3/4), and half (1/2) which meet flexibility and compatibility with a wide range of surgical instrument sets and medical equipment. Additionally, to meet specific needs, we offer various types of containers, providing particular solutions.





Colors



The containers offer six different color choices.



Ease of Identification

The containers feature color-coded lids, offering a quick and suitable method for categorization by color, which is easy to identify.





STANDARD MODEL CONTAINER SYSTEMS

Basic

12-14

Full Size (1/1), Quarter Size (3/4)
Half Size (1/2)

Long Body

15

Wide Body

16

Reusable Systems

18-21

Full Size (1/1), Quarter Size (3/4)
Half Size (1/2), Dental Size



Lid Perforated / Base Closed



Lid Safety Lid / Base Perforated



Lid & Base Perforated



Lid Safety Lid / Base Closed



STANDARD MODEL

Full Size (1/1)

Outside - Inside Dimensions (mm)

Code

585x280x90 - 550x260x85	C1-10-110
585x280x120 - 550x260x115	C1-10-135
585x280x140 - 550x260x135	C1-10-150
585x280x190 - 550x260x185	C1-10-200
585x280x250 - 550x260x245	C1-10-260



Base
Closed

585x280 (mm)

Code

●	C1-10-1
●	C1-10-2
●	C1-10-3
●	C1-10-4
●	C1-10-5
●	C1-10-6



Lid
Perforated

Outside - Inside Dimensions (mm)

Code

585x280x90 - 550x260x75	CP1-10-110
585x280x120 - 550x260x105	CP1-10-135
585x280x140 - 550x260x125	CP1-10-150
585x280x190 - 550x260x175	CP1-10-200
585x280x250 - 550x260x235	CP1-10-260



Base
Perforated *

585x280 (mm)

Code

●	CS1-10-1
●	CS1-10-2
●	CS1-10-3
●	CS1-10-4
●	CS1-10-5
●	CS1-10-6



Safety
Lid



STANDARD MODEL

Quarter Size (3/4)

Outside - Inside Dimensions (mm)

Code

465x280x90 - 430x260x85 C1-20-110

465x280x120 - 430x260x115 C1-20-135

465x280x140 - 430x260x135 C1-20-150

465x280x190 - 430x260x185 C1-20-180



Base
Closed

465x280 (mm)

Code

● C1-20-1

● C1-20-2

● C1-20-3

● C1-20-4

● C1-20-5

● C1-20-6



Lid
Perforated

Outside - Inside Dimensions (mm)

Code

465x280x90 - 430x260x75 CP1-20-110

465x280x120 - 430x260x105 CP1-20-135

465x280x140 - 430x260x125 CP1-20-150

465x280x190 - 430x260x175 CP1-20-180



Base
Perforated *

465x280 (mm)

Code

● CS1-20-1

● CS1-20-2

● CS1-20-3

● CS1-20-4

● CS1-20-5

● CS1-20-6



Safety
Lid

* Product code does not include the filters,
for visual purposes only.



STANDARD MODEL

Half Size (1/2)

Outside - Inside Dimensions (mm)

Code

295x280x90 - 260x260x85	C1-30-110
295x280x120 - 260x260x115	C1-30-135
295x280x140 - 260x260x135	C1-30-150
295x280x190 - 260x260x185	C1-30-200
295x280x250 - 260x260x245	C1-30-260



Base
Closed

295x280 (mm)

Code

●	C1-30-1
●	C1-30-2
●	C1-30-3
●	C1-30-4
●	C1-30-5
●	C1-30-6



Lid
Perforated

Outside - Inside Dimensions (mm)

Code

295x280x90 - 260x260x75	CPI-30-110
295x280x120 - 260x260x105	CPI-30-135
295x280x140 - 260x260x125	CPI-30-150
295x280x190 - 260x260x175	CPI-30-200
295x280x250 - 260x260x235	CPI-30-260



Base
Perforated *

295x280 (mm)

Code

●	CS1-30-1
●	CS1-30-2
●	CS1-30-3
●	CS1-30-4
●	CS1-30-5
●	CS1-30-6



Safety
Lid



STANDARD MODEL

Long Body

Outside - Inside Dimensions (mm)

Code

710x280x190 - 685x260x185

C1-40-190



Base
Closed

Outside - Inside Dimensions (mm)

Code

710x280x190 - 685x260x175

CP1-40-190



Base
Perforated *



Lid
Perforated

710x280 (mm)

Code

- ☐ C1-40-1
- ☒ C1-40-2
- ☐ C1-40-3
- ☐ C1-40-4
- ☐ C1-40-5
- ☐ C1-40-6

* Product code does not include the filters,
for visual purposes only.



STANDARD MODEL

Wide Body

Outside - Inside Dimensions (mm)

Code

600x400x110 - 565x380x105

C1-50-110

600x400x170 - 565x380x165

C1-50-170



Base
Closed

Outside - Inside Dimensions (mm)

Code

600x400x110 - 565x380x95

CP1-50-110

600x400x170 - 565x380x155

CP1-50-170



Base
Perforated *



Lid
Perforated

600x400 (mm)

Code

● C1-50-1

● C1-50-2

● C1-50-3

● C1-50-4


● C1-50-5

● C1-50-6



REUSABLE MODELS

Full Size (1/1)



Outside - Inside Dimensions (mm)

Code

585x280x90 - 550x260x85	C1-10-110
585x280x120 - 550x260x115	C1-10-135
585x280x140 - 550x260x135	C1-10-150
585x280x190 - 550x260x185	C1-10-200
585x280x250 - 550x260x245	C1-10-260



Base
Closed

Valve

585x280 mm

Code

- C2-10-1
- C2-10-2
- C2-10-3
- C2-10-4
- C2-10-5
- C2-10-6

Bio-Barrier

585x280 mm

Code

- C3-10-1
- C3-10-2
- C3-10-3
- C3-10-4
- C3-10-5
- C3-10-6

PTFE

585x280 mm

Code


- C4-10-1
- C4-10-2
- C4-10-3
- C4-10-4
- C4-10-5
- C4-10-6



Lid

REUSABLE MODELS

Quarter Size (3/4)



Outside - Inside Dimensions (mm)

Code

465x280x90 - 430x260x85	C1-20-110
465x280x120 - 430x260x115	C1-20-135
465x280x140 - 430x260x135	C1-20-150
465x280x190 - 430x260x185	C1-20-180



Base
Closed

Valve

465x280 (mm)

Code

- C2-20-1
- C2-20-2
- C2-20-3
- C2-20-4
- C2-20-5
- C2-20-6

Bio-Barrier

465x280 (mm)

Code

- C3-20-1
- C3-20-2
- C3-20-3
- C3-20-4
- C3-20-5
- C3-20-6

PTFE

465x280 (mm)

Code


- C4-20-1
- C4-20-2
- C4-20-3
- C4-20-4
- C4-20-5
- C4-20-6



Lid

REUSABLE MODELS

Half Size (1/2)



Outside - Inside Dimensions (mm)

Code

295x280x90 - 260x260x85	C1-30-110
295x280x120 - 260x260x115	C1-30-135
295x280x140 - 260x260x135	C1-30-150
295x280x190 - 260x260x185	C1-30-200
295x280x250 - 260x260x245	C1-30-260



Base
Closed

Valve

295x280 (mm)

Code

- C2-30-1
- C2-30-2
- C2-30-3
- C2-30-4
- C2-30-5
- C2-30-6

Bio Barrier

295x280 (mm)

Code

- C3-30-1
- C3-30-2
- C3-30-3
- C3-30-4
- C3-30-5
- C3-30-6

PTFE

295x280 (mm)

Code


- C4-30-1
- C4-30-2
- C4-30-3
- C4-30-4
- C4-30-5
- C4-30-6



Lid

REUSABLE MODELS

Dental Model



Outside - Inside Dimensions (mm)

Code

310x190x80 - 280x180x80	C1-60-80
310x190x100 - 280x180x100	C1-60-100
310x190x130 - 280x180x130	C1-60-130



Base
Closed

Valve

310x190 (mm)

Code

●	C2-60-1
●	C2-60-2
●	C2-60-3
●	C2-60-4
●	C2-60-5
●	C2-60-6

Bio-Barrier

310x190 (mm)

Code


●	C3-60-1
●	C3-60-2
●	C3-60-3
●	C3-60-4
●	C3-60-5
●	C3-60-6

PTFE


310x190 (mm)

Code

●	C4-60-1
●	C4-60-2
●	C4-60-3
●	C4-60-4
●	C4-60-5
●	C4-60-6



Lid



SMALL MODEL CONTAINER SYSTEMS

Flat & Mini Model

24-25

Dental & 1/2 Dental Model

26-27



Lid Perforated / Base Closed



Lid & Base Perforated

SMALL MODELS

Flat Model



Base
Closed

Outside - Inside Dimensions (mm)

Code

300x280x45 - 270x270x45

C1-90-45

300x280x75 - 270x270x75

C1-90-75



Base
Perforated *

Outside - Inside Dimensions (mm)

Code

300x280x75 - 270x270x65

CP1-90-75

300x280 (mm)

Code



C1-90-1



C1-90-2



C1-90-3



C1-90-4



C1-90-5




C1-90-6



Lid
Perforated

SMALL MODELS

Mini Model



Base
Closed

Outside - Inside Dimensions (mm)	Code
310x130x40 - 280x130x40	C1-70-40
310x130x70 - 280x130x70	C1-70-70
310x130x100 - 280x130x100	C1-70-100



Base
Perforated *

Outside - Inside Dimensions (mm)	Code
310x130x40 - 280x130x35	CP1-70-40
310x130x70 - 280x130x65	CP1-70-70
310x130x100 - 280x130x95	CP1-70-100




Lid
Perforated

310x140 (mm)	Code
	C1-70-1
	C1-70-2
	C1-70-3
	C1-70-4
	C1-70-5
	C1-70-6

* Product code does not include the filters,
for visual purposes only.

SMALL MODELS

Dental Model



Base
Closed



Base
Perforated *



Lid
Perforated

Outside - Inside Dimensions (mm)	Code
310x190x40 - 280x180x40	C1-60-40
310x190x60 - 280x180x60	C1-60-60
310x190x80 - 280x180x80	C1-60-80
310x190x100 - 280x180x100	C1-60-100
310x190x130 - 280x180x130	C1-60-130
310x190x40 - 280x180x35	CP1-60-40
310x190x60 - 280x180x55	CP1-60-60
310x190x80 - 280x180x75	CP1-60-80
310x190x100 - 280x180x95	CP1-60-100
310x190x130 - 280x180x125	CP1-60-130

310x190 (mm)	Code
	C1-60-1
	C1-60-2
	C1-60-3
	C1-60-4
	C1-60-5
	C1-60-6

SMALL MODELS

1/2 Dental Model



Base
Closed

Outside - Inside Dimensions (mm)	Code
210x155x50 - 180x145x50	C1-80-50
210x155x60 - 180x145x60	C1-80-60



Lid
Perforated

210x155 (mm)	Code
	C1-80-1
	C1-80-2
	C1-80-3
	C1-80-4
	C1-80-5
	C1-80-6

* Product code does not include the filters,
for visual purposes only.



SPECIFIC CONTAINERS

Implant & Scopy Model

30-31



Lid & Base Perforated

SPECIFIC

Implant Model



Types

Outside Dimensions (mm)
Inside Dimensions (mm)



Lid Perforated
Base Perforated

310x140x65
310x135x65

500x155x75
500x150x75



CP1-100-1

CP1-110-1



CP1-100-2

CP1-110-2



CP1-100-3

CP1-110-3



CP1-100-4

CP1-110-4



CP1-100-5

CP1-110-5



CP1-100-6

CP1-110-6



SPECIFIC

Scopy Model



Types

Outside Dimensions (mm)
Inside Dimensions (mm)



Lid Perforated
Base Perforated


250x60x50
250x55x45

450x70x70
450x65x65

CP1-120-1

CP1-130-1





PERFORATED TRAYS

Standard Model


34

Optic, Mini & Dental Model

35

Micro, Small, Large Model

36-37



Mesh Dimensions 4x4 mm
Range of Mesh 1 mm



Suitable for effortless adaptation of the components



Easy handling and carrying even when the lid is closed

PERFORATED TRAYS

Standard Model



Types

Outside Dimensions (mm)
Inside Dimensions (mm)



Polymer Foot
A10-10-1



Base (With Feet)



Lid

Half Size	Middle Size	DIN Size	Full Size	Long	Wide
B2-2525-30 242x252x40	B2-4025-30 405x252x40	B2-4825-30 484x252x40	B2-5425-30 539x252x40	B2-6625-30 655x252x40	B2-5333-30 530x330x40
B2-2525-50 242x252x60	B2-4025-50 405x252x60	B2-4825-50 484x252x60	B2-5425-50 539x252x60	B2-6625-50 655x252x60	B2-5333-50 530x330x60
B2-2525-70 242x252x80	B2-4025-70 405x252x80	B2-4825-70 484x252x80	B2-5425-70 539x252x80	B2-6625-70 655x252x80	B2-5333-70 530x330x80
B2-2525-100 242x252x110	B2-4025-100 405x252x110	B2-4825-100 484x252x110	B2-5425-100 539x252x110	B2-6625-100 655x252x110	B2-5333-100 530x330x110

B2-2525 246x256x15	B2-4025 409x256x15	B2-4825 488x256x15	B2-5425 543x256x15	B2-6625 659x256x15	B2-5333 534x334x15
-----------------------	-----------------------	-----------------------	-----------------------	-----------------------	-----------------------

PERFORATED TRAYS

Optic Model



Types

Outside Dimensions (mm)
Inside Dimensions (mm)

Single Scope



Base & Lid

B5-278-44
272x83x44

B5-458-44
452x83x44

Double Scope



Base & Lid

B5-2710-44
272x103x44

B5-4510-44
452x103x44



PERFORATED TRAYS

Mini & Dental Model

Mini Model



Base & Lid

B2-2712-45
270x123x45

B2-2712-65
270x123x65

Dental Model




Base & Lid

B2-2717-50
270x170x50

B2-2717-100
270x170x100


PERFORATED TRAYS WITH HINGES


Micro & Small Trays



Code	Outside Dimensions (mm)	Lock System	• Embossed lock system • 2x2mm perforated bottom and lid
B2-44-25	40x40x25	Embossed lock system	
B2-66-25	60x60x25	Embossed lock system	
B2-66-40	60x60x40	Embossed lock system	
B2-88-25	80x80x25	Embossed lock system	
B2-88-40	80x80x40	Embossed lock system	
B2-84-40	80x40x40	Embossed lock system	

Code	Outside Dimensions (mm)	Lock System	• Sliding flat lock system • 2x2mm perforated bottom and lid
B2-103-25	100x30x25	Push lock system	
B2-163-25	160x30x25	Push lock system	
B2-203-25	200x30x25	Push lock system	
B2-244-40	240x40x40	Push lock system	
B2-304-40	300x40x40	Push lock system	
B2-454-40	450x40x40	Push lock system	
B2-119-45	115x90x45	Sliding flat lock system	
B2-1111-25	115x115x25	Sliding flat lock system	
B2-1415-50	145x150x50	Sliding flat lock system	
B2-1713-35	170x135x35	Sliding flat lock system	
B2-2114-35	210x145x25	Sliding flat lock system	
B2-2114-45	210x145x45	Sliding flat lock system	
B2-2424-45	240x240x45	Sliding flat lock system	
B2-2717-35	275x175x35	Sliding flat lock system	
B2-2817-60	280x170x60	Sliding flat lock system	





Large Trays

Code	Outside Dimensions (mm)	Lock System	• Sliding flat lock system • 4x4mm perforated bottom and lid
B2-88-40-1	80x80x40	Push lock system	
B2-108-30	100x80x30	Sliding flat lock system	
B2-1311-40	130x115x40	Sliding flat lock system	
B2-1311-60	130x115x60	Sliding flat lock system	
B2-197-40	195x70x40	Sliding flat lock system	
B2-2213-40	225x130x40	Sliding flat lock system	
B2-2511-25	250x115x25	Sliding flat lock system	
B2-2511-40	250x115x40	Sliding flat lock system	
B2-2511-60	250x115x60	Sliding flat lock system	
B2-4811-50	480x115x50	Sliding flat lock system	
B2-2712-45-1	270x123x45	Sliding flat lock system	
B2-2712-65-1	270x123x65	Sliding flat lock system	
B2-2717-50-1	270x170x50	Sliding flat lock system	
B2-2717-100-1	270x170x100	Sliding flat lock system	

*Dimensions: width x length x height (with tolerances up to ± 5 mm)



HYBRID BASKETS

Standard Model


40-41

Specific Model

42



Base Plain Weave
Mesh Dimensions 4x4 mm
Wire Diameter 1 mm
Sides Perforated sheet
Perforation Diameter 3 mm



Base Welded Wire Mesh
Mesh Dimensions 6x6 mm
Wire Diameter 1.2 mm
Sides Perforated sheet
Perforation Diameter 3 mm



STANDARD

Plain Weave



Types

Outside Dimensions (mm)
Inside Dimensions (mm)



Base




Lid

Half Size	Middle Size	DIN Size	Full Size
B3-2525-50 245x255x50	B3-4025-50 405x255x50	B3-4825-50 485x255x50	B3-5425-50 540x255x50
B3-2525-60 245x255x60	B3-4025-60 405x255x60	B3-4825-60 485x255x60	B3-5425-60 540x255x60
B3-2525-70 245x255x70	B3-4025-70 405x255x70	B3-4825-70 485x255x70	B3-5425-70 540x255x70
B3-2525-100 245x255x100	B3-4025-100 405x255x100	B3-4825-100 485x255x100	B3-5425-100 540x255x100

Product Details

Base & Lid
Plain Weave



Side Plate
Perforated Sheet

STANDARD

Welded Wire Mesh



Types

Outside Dimensions (mm)
Inside Dimensions (mm)



Base

Half Size	Middle Size	DIN Size	Full Size
B3-2525-50-1 245x255x50	B3-4025-50-1 405x255x50	B3-4825-50-1 485x255x50	B3-5425-50-1 540x255x50
B3-2525-60-1 245x255x60	B3-4025-60-1 405x255x60	B3-4825-60-1 485x255x60	B3-5425-60-1 540x255x60
B3-2525-70-1 245x255x70	B3-4025-70-1 405x255x70	B3-4825-70-1 485x255x70	B3-5425-70-1 540x255x70
B3-2525-100-1 245x255x100	B3-4025-100-1 405x255x100	B3-4825-100-1 485x255x100	B3-5425-100-1 540x255x100




Lid

B3-2525-1 215x240	B3-4025-1 375x240	B3-4825-1 455x240	B3-5425-1 510x240
----------------------	----------------------	----------------------	----------------------

Product Details

Base & Lid
Welded Wire Mesh



Side Plate
Perforated Sheet

SPECIFIC

Plain Weave



Types

Outside Dimensions (mm)
Inside Dimensions (mm)




Base & Lid


Type 1	Type 2	Type 3
B3-1211-50 125x110x50	B3-2717-30 270x170x30	B3-3425-40 340x250x40
B3-1211-70 125x110x70	B3-2717-50 270x170x50	B3-4425-40 440x250x40
B3-197-40 195x70x40	B3-2717-70 270x170x70	B3-4425-80 440x250x80
B3-2212-40 220x125x40	B3-2717-100 270x170x100	B3-4632-80 460x320x80
B3-2612-30 265x120x30		
B3-2612-40 265x120x40		
B3-2612-50 265x120x50		
B3-2612-70 265x120x70		

Product Details

Base & Lid
Plain Weave



Side Plate
Perforated Sheet





WIRE BASKETS

Standard Model

46-50

Wide Body

51



Plain Weave

Mesh Dimensions 4x4 mm

Range of Mesh 1 mm



Welded Wire Mesh

Mesh Dimensions 6x6 mm

Range of Mesh 1.2 mm



Fine Mesh

Mesh Dimensions 2x2 mm

Range of Mesh 1 mm



STANDARD

Plain Weave



Types

Outside Dimensions (mm)
Inside Dimensions (mm)



Base



Lid

Half Size	Middle Size	DIN Size	Full Size
B1-2525-30 245x255x30	B1-4025-30 405x255x30	B1-4825-30 485x255x30	B1-5425-30 540x255x30
B1-2525-50 245x255x50	B1-4025-50 405x255x50	B1-4825-50 485x255x50	B1-5425-50 540x255x50
B1-2525-70 245x255x70	B1-4025-70 405x255x70	B1-4825-70 485x255x70	B1-5425-70 540x255x70
B1-2525-100 245x255x100	B1-4025-100 405x255x100	B1-4825-100 485x255x100	B1-5425-100 540x255x100

..... Product Details

Base & Lid
Plain Weave



STANDARD

Welded Wire Mesh



Types

Outside Dimensions (mm)
Inside Dimensions (mm)



Base




Lid

Half Size	Middle Size	DIN Size	Full Size
B1-2525-30-1 245x255x30	B1-4025-30-1 405x255x30	B1-4825-30-1 485x255x30	B1-5425-30-1 540x255x30
B1-2525-50-1 245x255x50	B1-4025-50-1 405x255x50	B1-4825-50-1 485x255x50	B1-5425-50-1 540x255x50
B1-2525-70-1 245x255x70	B1-4025-70-1 405x255x70	B1-4825-70-1 485x255x70	B1-5425-70-1 540x255x70
B1-2525-100-1 245x255x100	B1-4025-100-1 405x255x100	B1-4825-100-1 485x255x100	B1-5425-100-1 540x255x100

..... Product Details

Base & Lid
Welded Wire Mesh





STANDARD

Optic Wire Baskets

Types

Outside Dimensions (mm)
Inside Dimensions (mm)



Base & Lid

Single Scope

B5-298-55-S
290x80x55

B5-468-55-S
460x80x55

B5-678-55
670x80x55

Double Scope

B5-298-55-D
290x80x55

B5-468-55-D
460x80x55



Silicone Holder
(Type 1)



Base & Lid

Single Scope

B5-298-55-S1
290x80x55

B5-468-55-S1
460x80x55



Double Scope

B5-298-55-D1
290x80x55

B5-468-55-D1
460x80x55




Silicone Holder
(Type 2)



STANDARD

Mini & Dental Model



Types

Outside Dimensions (mm)
Inside Dimensions (mm)



Base

Mini Model
(Plain Weave)

B1-2611-43
265x120x45

Dental Model
(Plain Weave)

B1-2717-38
270x170x45



Fine Mesh Baskets



Base & Lid

B4-88-40
80x80x40

B4-1415-50
145x150x50


B4-2214-50
220x140x50

B4-2717-35
275x175x35

B4-2817-60
280x170x60



Base & Lid
Fine Wire Mesh



Wide Body Wire Baskets

Types

Outside Dimensions (mm)
Inside Dimensions (mm)




Base

Wide Body (Intermediate Crimp)

B1-4733-50
470x337x51

B1-5333-100
533x337x101




Surgical Laundry Baskets



Base

B6-6040-150
605x405x150

B6-6040-300
605x405x300



COMPONENTS

Container Components

54-55

Perforated Tray & Wire Basket Components

56-61

.....





Filter Retainer
(For Standard Model)

A6-10-1



Filter Retainer
(For Mini & Dental Containers)

A6-10-2



Filter Retainer
(For Mini Dental Container)

A6-10-3



Security Seal

A2-10-1
w/o indicator 100 pcs

A2-10-2
with indicator 100 pcs



Security Seal
(Service)

A2-20-1
w/o indicator 100 Pcs.

A2-20-2
with indicator 100 Pcs.

COMPONENTS

Container Components



Paper Filter
(For 1/1, 3/4, 1/2 Size)

A1-10-1
No Logo 100 Pcs.



Paper Filter
(For Mini & Dental Containers)

A1-12-1
No Logo 100 Pcs.



Paper Filter
(For Mini Dental Container)

A1-13-1
No Logo 100 Pcs.



PTFE Filter
(For 1/1, 3/4, 1/2 Size)

A1-20-1



PTFE Filter
(For Mini & Dental Containers)

A1-22-1



PTFE Filter
(For Mini Dental Container)

A1-23-1



Identification Labels

- A5-10-1
- A5-10-2
- A5-10-3
- A5-10-4
- A5-10-5
- A5-10-6



Indicators
(Paper)


- A3-10-1 Size Big - 100 Pcs.
- A3-10-2 Size Small - 100 Pcs

Silicone Mats

Code	Outside Dimensions
A7-2712	270x125 mm
A7-2223	225x235 mm
A7-3723	375x235 mm
A7-4523	450x235 mm (DIN)
A7-5023	500x235 mm




Silicone Mat
(Type 1 - Blue)



Code	Outside Dimensions
A7-2612-1	260x125 mm
A7-2223-1	225x235 mm
A7-3723-1	370x235 mm
A7-5023-1	500x235 mm




Silicone Mat
(Type 2 - Blue)



COMPONENTS

Holdings & Endoscopic Racks

B10-4823-160
Endoscopic Rack
485 x 235 x 160 mm




Silicone Holder
(For Optic Wire Baskets)

40x17x6,5 mm A8-10-1
70x17x6,5 mm A8-10-2
Complete A8-10-3




Silicone Holder
(For Optic Wire Baskets)

54x17x6,5 mm A8-10-4
68x17x6,5 mm A8-10-5
80,5x17x6,5 mm A8-10-6
Complete A8-10-7



Silicone Holder
(For Adjustable Racks)


226x17 mm A8-10-8
226x17 mm A8-10-9



230x220x40mm Complete B10-2322-40

1 pcs A8-10-8
2 pcs A8-10-9

*Dimensions: length x width x height



KERRISON Rack
(For 8 Kerrisons)

240x230x170 mm B10-2423-170
(Suitable for 1/2 Containers)



KERRISON Rack
(For 17 Kerrisons)

400x245x170 mm B10-4024-170
(Suitable for 3/4 Containers)



A8-10-10
17x25x40 mm

A8-10-12
17x100x40 mm

A8-10-14
17x225x40 mm



A8-10-11
17x25x40 mm

A8-10-13
17x100x40 mm

A8-10-15
17x225x40 mm



A8-10-16
17x25x40 mm

A8-10-18
17x100x40 mm

A8-10-20
17x225x40 mm



A8-10-17
17x25x40 mm

A8-10-19
17x100x40 mm


A8-10-21
17x225x40 mm



A8-10-22
17x25x40 mm

A8-10-24
17x100x40 mm

A8-10-26
17x225x40 mm



A8-10-23
17x25x40 mm

A8-10-25
17x100x40 mm

A8-10-27
17x225x40 mm



A8-10-28
12x25x3 mm

A8-10-29
12x100x3 mm

A8-10-30
12x225x3 mm



COMPONENTS

Silicone Holders



A8-10-31
17x25x37mm
Ø12mm



A8-10-32
17x25x37 mm
Ø14mm



A8-10-37
17x25x37 mm
Ø12mm



A8-10-38
17x25x37 mm
Ø14mm



A8-10-33
17x25x30 mm
Ø12mm



A8-10-34
17x25x30 mm
Ø14mm



A8-10-39
17x25x30 mm
Ø12mm



A8-10-40
17x25x30 mm
Ø14mm



A8-10-35
17x24x21 mm
Ø5mm




A8-10-36
17x24x21 mm
Ø10mm



A8-10-41
17x24x21 mm
Ø5mm



A8-10-42
17x24x21 mm
Ø10mm





A8-10-43
17x225x21 mm
Ø10mm - Ø5mm

It can be used with endoscopic rack "B10-4823-160" if desired



*Dimensions: width x length x height

Hanging Trays




Base

Code	Outside Dimensions
B2-2570-30	254x70x30 mm
B2-2517-30	254x120x30 mm
B2-2512-30	254x170x30 mm
B2-2524-30	254x240x30 mm

Compatible with standard tray and basket series.




Individually usable for sterilization requirements.




COMPONENTS

Separators

Non-Sliding



Sliding



Code

Outside Dimensions

A9-10-1	230x240x45mm
A9-10-2	390x240x45mm
A9-10-3	470x240x45mm
A9-10-4	530x240x45mm

Code

Outside Dimensions

A9-10-5	230x240x45mm
A9-10-6	390x240x45mm
A9-10-7	470x240x45mm
A9-10-8	530x240x45mm

Type	Packing		Steam sterilization procedure acc. To DIN EN 868-8:2009-09				
	Description	Symbolic Delivery and Packaging System	Gravity Procedure	Pre-Vacuum Procedure	Fractional Vacuum Procedure	Fractional Flow Procedure	Steam Injection Procedure
With filters (StB-F)	Textiles in sterilization container with filter in lid and bottom acc. To DIN EN 868-8:2009-09		+	+	+	+	+
			0	0	+	+	+
	Textiles in sterilization container with filter in lid, bottom not perforated to DIN EN 868-8:2009-09		-	-	+	+	+
			-	-	+	+	0
	Textiles in sterilization container with filter in bottom, lid not perforated to DIN EN 868-8:2009-09		0	+	+	+	+
			-	0	+	+	+
	Instruments in wire basket acc. To DIN 58952 part 3L in sterilization container with filter in lid and bottom acc. To DIN EN 868-8:2009-09		+	+	+	+	+
	Instruments in wire basket acc. To DIN 58952 part 3L in sterilization container with filter in lid and bottom not perforated acc. To DIN EN 868-8:2009-09		-	+	+	+	+
	Instruments in wire basket acc. To DIN 58952 part 3L in sterilization container with filter in bottom acc. To DIN EN 868-8:2009-09		0	+	+	+	+
			0	+	+	+	+
With valve (StB-V)	Textiles in sterilization container with valves in lid and bottom acc. To DIN EN 868-8:2009-09		-	0	+	+	0
			-	-	+	+	0
	Textiles in sterilization container with valves in lid, no valves on the bottom		-	-	+	+	-
			-	-	+	+	-
	Instruments in wire basket in sterilization container with valve in lid and bottom		-	+	+	+	+
	Instruments in wire basket in sterilization container with valve in lid no valves on the bottom		-	+	+	+	+
			-	+	+	+	+

INDEX

Numerical

Code	Page	Code	Page	Code	Page	Code	Page	Code	Page
A1-10-2.....	55	A8-10-18.....	58	B1-4025.....	46	B2-1311-40.....	37	B2-5425-50.....	34
A1-12-2.....	55	A8-10-19.....	58	B1-4025-1.....	47	B2-1311-60.....	37	B2-5425-70.....	34
A1-13-2.....	55	A8-10-20.....	58	B1-4025-30.....	46	B2-1415-50.....	36	B2-5425-100.....	34
A1-20-1.....	55	A8-10-21.....	58	B1-4025-30-1.....	47	B2-1713-35.....	36	B2-6625.....	34
A1-22-1.....	55	A8-10-22.....	58	B1-4025-50.....	46	B2-2114-35.....	36	B2-6625-30.....	34
A1-23-1.....	55	A8-10-23.....	58	B1-4025-50-1.....	47	B2-2114-45.....	36	B2-6625-50.....	34
A2-10-1.....	54	A8-10-24.....	58	B1-4025-100.....	46	B2-2213-40.....	37	B2-6625-70.....	34
A2-20-1.....	54	A8-10-25.....	58	B1-4025-100-1.....	47	B2-2424-45.....	36	B2-6625-100.....	34
A3-10-1.....	55	A8-10-26.....	58	B1-4733-50.....	51	B2-2511-25.....	37	B3-197-40.....	42
A3-10-2.....	55	A8-10-27.....	58	B1-4825.....	46	B2-2511-40.....	37	B3-1211-50.....	42
A5-10-1.....	55	A8-10-28.....	58	B1-4825-1.....	47	B2-2511-60.....	37	B3-1211-70.....	42
A5-10-2.....	55	A8-10-29.....	58	B1-4825-30.....	46	B2-2512-30.....	60	B3-2212-40.....	42
A5-10-3.....	55	A8-10-30.....	58	B1-4825-30-1.....	47	B2-2517-30.....	60	B3-2525.....	40
A5-10-4.....	55	A8-10-31.....	59	B1-4825-50.....	46	B2-2524-30.....	60	B3-2525-1.....	41
A5-10-5.....	55	A8-10-32.....	59	B1-4825-50-1.....	47	B2-2525.....	34	B3-2525-50.....	40
A5-10-6.....	55	A8-10-33.....	59	B1-4825-70.....	46	B2-2525-30.....	34	B3-2525-50-1.....	41
A6-10-1.....	54	A8-10-34.....	59	B1-4825-70-1.....	47	B2-2525-50.....	34	B3-2525-60.....	40
A6-10-2.....	54	A8-10-35.....	59	B1-4825-100.....	46	B2-2525-70.....	34	B3-2525-60-1.....	41
A6-10-3.....	54	A8-10-36.....	59	B1-4825-100-1.....	47	B2-2525-100.....	34	B3-2525-70.....	40
A7-2223.....	56	A8-10-37.....	59	B1-5333-100.....	51	B2-2570-30.....	60	B3-2525-70-1.....	41
A7-2223-1.....	56	A8-10-38.....	59	B1-5425.....	46	B2-2712-45.....	35	B3-2525-100.....	40
A7-2612-1.....	56	A8-10-39.....	59	B1-5425-1.....	47	B2-2712-45-1.....	37	B3-2525-100-1.....	41
A7-2712.....	56	A8-10-40.....	59	B1-5425-30.....	46	B2-2712-65.....	35	B3-2612-30.....	42
A7-3723.....	56	A8-10-41.....	59	B1-5425-30-1.....	47	B2-2712-65-1.....	37	B3-2612-40.....	42
A7-3723-1.....	56	A8-10-42.....	59	B1-5425-50.....	46	B2-2717-35.....	36	B3-2612-50.....	42
A7-4523.....	56	A8-10-43.....	59	B1-5425-50-1.....	47	B2-2717-50.....	35	B3-2612-70.....	42
A7-5023.....	56	A9-10-1.....	61	B1-5425-70.....	46	B2-2717-50-1.....	37	B3-2717-30.....	42
A7-5023-1.....	56	A9-10-2.....	61	B1-5425-70-1.....	47	B2-2717-100.....	35	B3-2717-50.....	42
A8-10-1.....	57	A9-10-3.....	61	B1-5425-100.....	46	B2-2717-100-1.....	37	B3-2717-100.....	42
A8-10-2.....	57	A9-10-4.....	61	B1-5425-100-1.....	47	B2-2817-60.....	36	B3-3425-40.....	42
A8-10-3.....	57	A9-10-5.....	61	B2-44-25.....	36	B2-4025.....	34	B3-4025.....	40
A8-10-4.....	57	A9-10-6.....	61	B2-66-25.....	36	B2-4025-30.....	34	B3-4025-1.....	41
A8-10-5.....	57	A9-10-7.....	61	B2-66-40.....	36	B2-4025-50.....	34	B3-4025-50.....	40
A8-10-6.....	57	A9-10-8.....	61	B2-84-40.....	36	B2-4025-100.....	34	B3-4025-50-1.....	41
A8-10-7.....	57	A10-10-1.....	34	B2-88-25.....	36	B2-4811-50.....	37	B3-4025-60.....	40
A8-10-8.....	57	B1-2525.....	46	B2-88-40.....	36	B2-4825.....	34	B3-4025-60-1.....	41
A8-10-8.....	57	B1-2525-1.....	47	B2-88-40-1.....	37	B2-4825-30.....	34	B3-4025-100.....	40
A8-10-9.....	57	B1-2525-30.....	46	B2-103-25.....	36	B2-4825-50.....	34	B3-4025-100-1.....	41
A8-10-9.....	57	B1-2525-30-1.....	47	B2-108-30.....	37	B2-4825-70.....	34	B3-4425-40.....	42
A8-10-10.....	58	B1-2525-50.....	46	B2-119-45.....	36	B2-4825-100.....	34	B3-4425-80.....	42
A8-10-11.....	58	B1-2525-50-1.....	47	B2-163-25.....	36	B2-5333.....	34	B3-4632-80.....	42
A8-10-12.....	58	B1-2525-70.....	46	B2-197-40.....	37	B2-5333-30.....	34	B3-4825.....	40
A8-10-13.....	58	B1-2525-70-1.....	47	B2-203-25.....	36	B2-5333-50.....	34	B3-4825-1.....	41
A8-10-14.....	58	B1-2525-100.....	46	B2-244-40.....	36	B2-5333-70.....	34	B3-4825-50.....	40
A8-10-15.....	58	B1-2525-100-1.....	47	B2-304-40.....	36	B2-5333-100.....	34	B3-4825-50-1.....	41
A8-10-16.....	58	B1-2611-43.....	50	B2-454-40.....	36	B2-5425.....	34	B3-4825-60.....	40
A8-10-17.....	58	B1-2717-38.....	50	B2-1111-25.....	36	B2-5425-30.....	34	B3-4825-60-1.....	41

Code	Page	Code	Page	Code	Page	Code	Page	Code	Page
B3-4825-70	40	C1-20-1	13	C1-60-80	26	C3-10-1	18	CPI-20-135	13
B3-4825-70-1	41	C1-20-2	13	C1-60-100	21	C3-10-2	18	CPI-20-150	13
B3-4825-100	40	C1-20-3	13	C1-60-100	26	C3-10-3	18	CPI-20-180	13
B3-4825-100-1	41	C1-20-4	13	C1-60-130	21	C3-10-4	18	CPI-30-110	14
B3-5425	40	C1-20-5	13	C1-60-130	26	C3-10-5	18	CPI-30-135	14
B3-5425-1	41	C1-20-6	13	C1-70-1	25	C3-10-6	18	CPI-30-150	14
B3-5425-50	40	C1-20-110	13	C1-70-2	25	C3-20-1	19	CPI-30-200	14
B3-5425-50-1	41	C1-20-110	19	C1-70-3	25	C3-20-2	19	CPI-30-260	14
B3-5425-60	40	C1-20-135	13	C1-70-4	25	C3-20-3	19	CPI-40-190	15
B3-5425-60-1	41	C1-20-135	19	C1-70-5	25	C3-20-4	19	CPI-50-110	16
B3-5425-70	40	C1-20-150	13	C1-70-6	25	C3-20-5	19	CPI-50-170	16
B3-5425-70-1	41	C1-20-150	19	C1-70-40	25	C3-20-6	19	CPI-60-40	26
B3-5425-100	40	C1-20-180	13	C1-70-70	25	C3-30-1	20	CPI-60-60	26
B3-5425-100-1	41	C1-20-180	19	C1-70-100	25	C3-30-2	20	CPI-60-80	26
B4-88-40	50	C1-30-1	14	C1-80-1	27	C3-30-3	20	CPI-60-100	26
B4-1415-50	50	C1-30-2	14	C1-80-2	27	C3-30-4	20	CPI-60-130	26
B4-2214-50	50	C1-30-3	14	C1-80-3	27	C3-30-5	20	CPI-70-40	25
B4-2717-35	50	C1-30-4	14	C1-80-4	27	C3-30-6	20	CPI-70-70	25
B4-2817-60	50	C1-30-5	14	C1-80-5	27	C3-60-1	21	CPI-70-100	25
B5-278-44	35	C1-30-6	14	C1-80-6	27	C3-60-2	21	CPI-90-75	24
B5-298-55-D	48	C1-30-110	14	C1-80-50	27	C3-60-3	21	CPI-100-1	30
B5-298-55-D1	48	C1-30-110	20	C1-80-60	27	C3-60-4	21	CPI-100-2	30
B5-298-55-S	48	C1-30-135	14	C1-90-1	24	C3-60-5	21	CPI-100-3	30
B5-298-55-S1	48	C1-30-135	20	C1-90-2	24	C3-60-6	21	CPI-100-4	30
B5-458-44	35	C1-30-150	14	C1-90-3	24	C4-10-1	18	CPI-100-5	30
B5-468-55-D	48	C1-30-150	20	C1-90-4	24	C4-10-2	18	CPI-100-6	30
B5-468-55-D1	48	C1-30-200	14	C1-90-5	24	C4-10-3	18	CPI-110-1	30
B5-468-55-S	48	C1-30-200	20	C1-90-6	24	C4-10-4	18	CPI-110-2	30
B5-468-55-S1	48	C1-30-260	14	C1-90-45	24	C4-10-5	18	CPI-110-4	30
B5-678-55	48	C1-30-260	20	C1-90-75	24	C4-10-6	18	CPI-110-5	30
B5-2710-44	35	C1-40-1	15	C2-10-1	18	C4-20-1	19	CPI-110-6	30
B5-4510-44	35	C1-40-2	15	C2-10-2	18	C4-20-2	19	CPI-120-1	31
B6-6040-150	51	C1-40-3	15	C2-10-3	18	C4-20-3	19	CPI-130-1	31
B6-6040-300	51	C1-40-4	15	C2-10-4	18	C4-20-4	19	CSI-10-1	12
B10-2322-40	57	C1-40-5	15	C2-10-5	18	C4-20-5	19	CSI-10-2	12
B10-2423-170	57	C1-40-6	15	C2-10-6	18	C4-20-6	19	CSI-10-3	12
B10-4024-170	57	C1-40-190	15	C2-20-1	19	C4-30-1	20	CSI-10-4	12
B10-4823-160	57	C1-50-1	16	C2-20-2	19	C4-30-2	20	CSI-10-5	12
C1-10-1	12	C1-50-2	16	C2-20-3	19	C4-30-3	20	CSI-10-6	12
C1-10-2	12	C1-50-3	16	C2-20-4	19	C4-30-4	20	CSI-20-1	13
C1-10-3	12	C1-50-4	16	C2-20-5	19	C4-30-5	20	CSI-20-2	13
C1-10-4	12	C1-50-5	16	C2-20-6	19	C4-30-6	20	CSI-20-3	13
C1-10-5	12	C1-50-6	16	C2-30-1	20	C4-60-1	21	CSI-20-4	13
C1-10-6	12	C1-50-110	16	C2-30-2	20	C4-60-2	21	CSI-20-5	13
C1-10-110	12	C1-50-170	16	C2-30-3	20	C4-60-3	21	CSI-20-6	13
C1-10-110	18	C1-60-1	26	C2-30-4	20	C4-60-4	21	CSI-30-1	14
C1-10-135	12	C1-60-2	26	C2-30-5	20	C4-60-5	21	CSI-30-2	14
C1-10-135	18	C1-60-3	26	C2-30-6	20	C4-60-6	21	CSI-30-3	14
C1-10-150	12	C1-60-4	26	C2-60-1	21	CPI-10-110	12	CSI-30-4	14
C1-10-150	18	C1-60-5	26	C2-60-2	21	CPI-10-135	12	CSI-30-5	14
C1-10-200	12	C1-60-6	26	C2-60-3	21	CPI-10-150	12	CSI-30-6	14
C1-10-200	18	C1-60-40	26	C2-60-4	21	CPI-10-200	12		
C1-10-260	12	C1-60-60	26	C2-60-5	21	CPI-10-260	12		
C1-10-260	18	C1-60-80	21	C2-60-6	21	CPI-20-110	13		