

# ELECTROSURGERY

## ELECTROCAUTERIES

FIAB sterile cautery pens are used in haemostasis, granting a fast and precise cauterization of small blood vessels. They are extremely easy to handle, thanks to an internal long life battery; several models are available, different from each other for working temperatures (1200°C, 800°C and 600°C), shape, thickness and length of the cautery tip.

The available range allows to choose a cautery pen that perfectly matches the clinical demands of the procedure.



DISPOSABLE CAUTERIES							
REF	F7244	F7266	F7234	F7277	F7288	F7255	F7255UF
TIP SHAPE	FINE	LARGE	FINE	FINE	THICK	FINE	ULTRA FINE
TIP LENGTH (mm )	28	28	75	125	28	28	28
TEMPERATURE AT THE ACTIVATION (°C / °F)	1200 / 2200				800 / 1400	600 / 1100	
POWER SUPPLY (V)	3,0				3,0	1,5	
DIAMETER (MM)	16				16	16	
WEIGHT (GR)	65	65	78	78	65	45	45
SERIES	HTC				MTC	LTC	

HIGH TEMPERATURE  
CAUTERY

MED. TEMPERATURE  
CAUTERY

LOW TEMPERATURE  
CAUTERY



LARGE TIP



FINE TIP



THICK TIP



ULTRA FINE TIP

# ELECTROSURGERY

F7297



F7299

REUSABLE CAUTERIES						
REF	F7297	F7298F	F7299	F7298F	F7298L	F7298FG
TIP SHAPE	--	LARGE	--	FINE	FINE	FINE
TIP LENGTH (mm )	12	3	18	3	3	12
TEMPERATURE AT THE ACTIVATION (°C / °F)	600/1100	--	1200 / 2200	--	--	--
POWER SUPPLY (V)	1,5	--	3,0	--	--	--
DIAMETER (MM)	16	--	16	--	--	--
WEIGHT (GR)	48	2	65	2	2	8
SERIES	LTC		HTC			

LOW TEMPERATURE  
CAUTERY AND TIP

HIGH TEMPERATURE  
CAUTERY AND TIP

F7298F

F7298L

F7298FG



## NERVE LOCATOR

### Neuropacer

**F1745** Disposable, sterile, battery operating, nerve locator complete with active electrode and grounding needle

